Spectroscopie Raman

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RAMAN DATA AND ANALYSIS

Raman Spectroscopy for Analysis and Monitoring

The Raman scattering technique is a vibrational molecular spectroscopy which derives from an inelastic light scattering process. With Raman spectroscopy, a laser photon is scattered by a sample molecule and loses (or gains) energy during the process. The amount of energy lost is seen as a change in energy (wavelength) of the irradiating photon. This energy loss is characteristic for a particular bond in the molecule. Raman can best be thought of as producing a precise spectral fingerprint, unique to a molecule or indeed and individual molecular structure. In this respect it is similar to the more commonly found FT-IR spectroscopy. However, unlike FT-IR, there are a distinct number of advantages when using Raman.

- Raman can be used to analyse aqueous solutions since it does not suffer from the large water absorption effects found with FT techniques.*
- The intensity of spectral features in solution is directly proportional to the concentration of the particular species*
- Raman spectra are generally robust to temperature changes*
- Raman requires little or no sample preparation. It does not need the use of Nujol, or KBr matrices and is largely unaffected by sample cell materials such as glass.*
- The use of a Raman microscope such as the LabRAM provides very high level of spatial resolution and depth discrimination, not found with the FT methods of analysis*

These advantages and its highly specific nature, mean that Raman has become a very powerful tool for analysis and chemical monitoring. Depending upon instrumentation, it is a technique which can be used for the analysis of solids, liquids and solutions and can even provide information on physical characteristics such as crystalline phase and orientation, polymorphic forms, and intrinsic stress.

Functional Group/ Vibration	Region	Raman	InfraRed
Lattice vibrations in crystals, LA modes	10 - 200 cm ⁻¹	strong	strong
$\delta(\text{CC})$ aliphatic chains	250 - 400 cm ⁻¹	strong	weak
$\nu(\text{Se-Se})$	290 - 330 cm ⁻¹	strong	weak
$\nu(\text{S-S})$	430 - 550 cm ⁻¹	strong	weak
$\nu(\text{Si-O-Si})$	450 - 550 cm ⁻¹	strong	weak
$\nu(\text{Xmetal-O})$	150-450 cm ⁻¹	strong	med-weak
$\nu(\text{C-I})$	480 - 660 cm ⁻¹	strong	strong
$\nu(\text{C-Br})$	500 - 700 cm ⁻¹	strong	strong
$\nu(\text{C-Cl})$	550 - 800 cm ⁻¹	strong	strong
$\nu(\text{C-S})$ aliphatic	630 - 790 cm ⁻¹	strong	medium
$\nu(\text{C-S})$ aromatic	1080 - 1100 cm ⁻¹	strong	medium
$\nu(\text{O-O})$	845 - 900 cm ⁻¹	strong	weak
$\nu(\text{C-O-C})$	800 - 970 cm ⁻¹	medium	weak
$\nu(\text{C-O-C})$ asym	1060 - 1150 cm ⁻¹	weak	strong

			Raman Bands
$\nu(\text{CC})$ alicyclic, aliphatic chain vibrations	600 - 1300 cm^{-1}	medium	Medium
$\nu(\text{C=S})$	1000 - 1250 cm^{-1}	strong	weak
$\nu(\text{CC})$ aromatic ring chain vibrations	*1580, 1600 cm^{-1}	strong	medium
	**1450, 1500 cm^{-1}	medium	medium
	*1000 cm^{-1}	strong/medium	weak
$\delta(\text{CH}_3)$	1380 cm^{-1}	medium	strong
$\delta(\text{CH}_2)$			
$\delta(\text{CH}_3)$ asym	1400 - 1470 cm^{-1}	medium	medium
$\delta(\text{CH}_2)$	1400 - 1470 cm^{-1}	medium	medium
$\delta(\text{CH}_3)$ asym	1400 - 1470 cm^{-1}	medium	medium
$\nu(\text{C-(NO}_2))$	1340 - 1380 cm^{-1}	strong	medium
$\nu(\text{C-(NO}_2))$ asym	1530 - 1590 cm^{-1}	medium	strong
$\nu(\text{N=N})$ aromatic	1410 - 1440 cm^{-1}	medium	-
$\nu(\text{N=N})$ aliphatic	1550 - 1580 cm^{-1}	medium	-
$\delta(\text{H}_2\text{O})$	~1640 cm^{-1}	weak broad	strong
$\nu(\text{C}=\text{N})$	1610 - 1680 cm^{-1}	strong	medium
$\nu(\text{C}=\text{C})$	1500 - 1900 cm^{-1}	strong	weak
$\nu(\text{C}=\text{O})$	1680 - 1820 cm^{-1}	medium	strong
$\nu(\text{C}=\text{C})$	2100 - 2250 cm^{-1}	strong	weak
$\nu(\text{C}=\text{N})$	2220 - 2255 cm^{-1}	medium	strong
$\nu(\text{-S-H})$	2550 - 2600 cm^{-1}	strong	weak
$\nu(\text{C-H})$	2800 - 3000 cm^{-1}	strong	strong
$\nu(\text{=C-H})$	3000 - 3100 cm^{-1}	strong	medium
$\nu(\text{≡C-H})$	3300 cm^{-1}	weak	strong
$\nu(\text{N-H})$	3300 - 3500 cm^{-1}	medium	medium
$\nu(\text{O-H})$	3100 - 3650 cm^{-1}	weak	strong

HORIBA JOBIN YVON

(All HORIBA Jobin Yvon companies were formerly known as Jobin Yvon)

France : HORIBA Jobin Yvon S.A.S., 231 rue de Lille, 59650 Villeneuve d'Ascq. Tel : +33 (0)3 20 59 18 00.
USA : HORIBA Jobin Yvon Inc., 3880 Park Avenue, Edison, NJ 08820-3012. Tel : +1-732-494-8660,
Fax : +1-732-549-2571. Email : raman@jobinyvon.com www.jobinyvon.com
Japan : HORIBA Ltd., JY Optical Sales Dept., 1-7-8 Higashi-kanda, Chiyoda-ku, Tokyo 101-0031.
Tel : +81 (0)3 3861 8231. Fax : +81 (0)3 3861 8259. Email: raman@horiba.com
Germany: +49 (0) 6251 84 75-0 **Italy:** +39 02 57603050 **UK:** +44 (0)20 8204 8142
China: +86 (0) 10 6849 2216

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Determination of Unsaturation in Food Oils and Fats

RAMANRXN SYSTEMS™ Application Note

Number 500

Key Issues

- Straightforward determination of total unsaturation and *cis* isomer content
- High water content, preventing use of FTIR
- Rapid on-line process control

Introduction

Knowing the degree of unsaturation and relative amounts of *cis* and *trans* isomerization in edible oils and fats is important in today's foods industry. This is necessary for both process and quality control, and is required nutritional information in food labeling.

Total oil and fat unsaturation is presented as the Iodine Value (IV). The IV is currently determined by traditional titration methods. FTIR is the standard technique to determine *trans* isomer content. GC is used to estimate lipid fatty acid composition. However, IR absorption by water and chromatographic time constraints, along with stringent sampling requirements, make these techniques unsuitable for on-line process control.

Raman spectroscopy has been proposed as a technique suited for the rapid on-line determination of both *cis/trans* isomer content and total unsaturation in oils and fats.¹⁻³ An off-line demonstration of this analysis has appeared using FT-Raman.³ We present evidence that this analysis can be conducted with dispersive Raman instrumentation designed for process control.

Experimental

The data was acquired using a fully integrated, remote Kaiser Optical Systems, Inc. HoloProbe™ analyzer with 10 mW of 532-nm incident laser light. A HoloPlex™ grating was used to acquire the full Raman spectrum in one exposure, allowing quantitative analysis across the entire measured range. Integration times ranged from 60 to 300 s, depending on the sample.

Determination of Total Unsaturation

Figure 1 shows the Raman spectra of butter, margarine, vegetable oil and vegetable shortening

between 1150 and 1775 cm⁻¹. Total unsaturation can be determined quantitatively by measuring the ratio of the C=C stretch (ν C=C) centered at 1661 cm⁻¹ to the CH₂ scissoring deformation (δ CH₂) at 1444 cm⁻¹.¹⁻³ When comparisons are made between the relative intensities of ν C=C and δ CH₂ in Figure 1, it is clear that butter has the highest percentage of saturated fat. Conversely, vegetable oil has the least saturation. Vegetable shortening and margarine fall between these two extremes. This matches what we know qualitatively. Vegetable oils and shortenings are nutritionally preferable to animal shortenings because they contain less saturated fat.


Figure 1. Background-corrected Raman spectra of commercial food oils and fats. The bands at 1165 and 1525 cm⁻¹ of the butter and margarine spectra are attributed to the coloring agent β-carotene.

Cis Isomer Content

Cis isomer content can also be determined from the Raman spectrum. The band at 1272 cm⁻¹ is attributed to in-plane =C–H deformation in an unconjugated *cis* double bond. The intensity of this band ratioed to the in-phase methylene twisting

RAMANRXN SYSTEMS
REACTION ANALYSIS, MONITORING AND CONTROL

vibration at 1306 cm^{-1} provides a direct *cis* isomer measurement. This is illustrated by comparing the Raman spectrum of vegetable oil with that of partially hydrogenated vegetable shortening. Hydrogenation results in a preferential decrease in the *cis/trans* isomer ratio.³ In Figure 1, unhydrogenated vegetable oil contains significantly more *cis* isomer character than partially hydrogenated vegetable shortening.

Background Correction in Highly Fluorescent Samples

Figure 2 shows the background-corrected, full Raman spectra of vegetable and peanut oils. The ratio of background level to the δCH_2 vibration intensity appears in parentheses by each spectrum. The fluorescence observed from commercial samples has been attributed to assorted coloring agents.¹⁻³ Figure 2 illustrates that, even in the presence of significant background, Raman spectra from these sample classes are readily observable.

In those instances where fluorescence is not overwhelming, dispersive Raman with visible excitation offers distinct advantages over FT-Raman.

Sensitivity is enhanced due to both the v^4 Raman intensity dependence and the more efficient CCD detector. Alternatively, fluorescence can often be significantly reduced or eliminated by using a longer incident wavelength. A **RAMANRXN1™** analyzer

equipped with a 785-nm NIR laser, is well-suited for this purpose.

Conclusion

The high-throughput **RAMANRXN1** analyzer permits the use of low laser powers and short exposure times for remote Raman measurements. Its small size, rugged construction, and minimal utility needs make it well suited to applications outside the laboratory environment. This technique is well suited to the online determination of total unsaturation and *cis* isomer level in commercial oil and fat samples.

*Current model is **RAMANRXN1** analyzer

References:

1. Sadeghi-Jorabchi, H.; Hendra, P. J.; Wilson, R. H.; Belton, P. S.; Journal of the American Oil Chemistry Society, Vol. 67, 1990, 483.
2. Goral, J.; Zichy, V.; Spectrochimica Acta, Vol. 46A, 1990, 253.
3. Sadeghi-Jorabchi, H.; Wilson, R. H.; Belton, P. S.; Spectrochimica Acta, Vol. 47A, 1991, 1449.
4. Lewis, I.R. Process Raman Spectroscopy. In *Handbook of Raman Spectroscopy*; Lewis, I.R., Edwards, H.G.M., Eds. Marcel Dekker: New York, 2001; pp. 919-973.


Figure 2. Background-corrected spectra of vegetable and peanut oils.

RamanRxn Systems, RamanRxn1, and HoloPlex are trademarks of Kaiser Optical Systems, Inc. All rights reserved.
Copyright © 2001 Kaiser Optical Systems, Inc. All rights reserved.

WWW.KOSI.COM

371 Parkland Plaza
Ann Arbor, MI 48103
phone (734)665-8083 / fax (734)665-8199